DE

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PARASORB® Cone, PARASORB® Fleece HD

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PARASORB® Cone, PARASORB® Fleece HD

EN Please read instructions carefully prior to use.

APPLICATION FORM

The sterile absorbable Collagens PARASORB® Cone and Fleece HD are used as haemostatic agents after tooth extractions and oral surgery. The application of PARASORB® Cone and Fleece HD into the defect cavity supports the stabilization of the produced blood clot whilst helping to control bleeding at the application site.

COMPOSITION

1 PARASORB® Cone contains 22.4 mg of native collagen fibrils of equine origin. This is the same concentration as PARASORB® Fleece HD.

- 1 cm² PARASORB® Fleece HD contains 5.6 mg of native collagen fibrils of equine origin;
- 1 PARASORB® Fleece HD 1.8 x 3.6 cm contains 36.2 mg native collagen fibrils of equine origin.
- 1 PARASORB® Fleece HD 3 x 7 cm contains 117.6 mg native collagen fibrils of equine origin.

INTENDED USE

The PARASORB® Cone and Fleece HD are used as haemostatic agents after tooth extractions and oral surgery. PARASORB® Cone and Fleece HD can be combined with fibrin glue in patients with an increased risk of bleedings.

MODE OF APPLICATION

- After removal from the sterile packaging PARASORB® Cone and Fleece HD can be cut to the desired shape and size with sterile instruments.

NOTES FOR POST-OPERATIVE CARE FOLLOWING DENTAL SURGERY

To be avoided:

- Rinse your mouth out on the day of surgery (may dislodge clot).
- Drink alcohol or hot drinks on the day of surgery (may increase blood pressure or induce hyperaemia).
- Smoke for at least 24 hours (causes vasoconstriction).
- Carry out physical activity on the day of surgery (increasing blood pressure).

Recommended:

- Sip lukewarm drinks carefully (helps avoid clot removal by washing out socket).
- Bite on gauze or clean, damp handkerchief for 10 minutes if area bleeds, preferably sitting upright in a chair.
- Rinse mouth gently with hot salt mouthwash or Chlorhexidine (0.2%) mouthwash after 24 hours. Continue regular mouthwash use over the following week, especially after meals.
- Clean teeth as normal, avoiding area of surgery for 24 hours.

STERILITY

PARASORB® Cone and Fleece HD may not be resterilised and opened single packs must be discarded. Do not re-use PARASORB® Cone and Fleece HD in whole or in part, as it may compromise sterility and/or the performance of the PARASORB® Cone and Fleece HD.

The PARASORB® Cone and Fleece HD are used in dry state with light pressure to achieve better adhesion and haemostasis.

The PARASORB® Cone and Fleece HD is inserted into the extraction socket or the cavity as required.

Once inserted into the desired cavity PARASORB® Cone and Fleece HD can be sutured over and/or held in place with fibrin glue as required.

Due to the affinity of collagen with blood covered surfaces, dry instruments and gloves should be used when administering the PARASORB® Cone and Fleece HD.

Due to the absorbability of the PARASORB® Cone and Fleece HD they do not have to be removed. PARASORB® Cone and Fleece HD will be absorbed within 3 – 4 months; absorption time may differ as it is dependent upon multiple patient specific factors including but not limited to site of implantation, comorbidities and blood perfusion.

PARASORB® Cone and Fleece HD may only be used by qualified medical professionals.

Medical treatment of severe bleeding disorders is mandatory prior to the beginning of the procedure.

Unless otherwise prescribed, the PARASORB® Cone and Fleece HD can be applied directly into the extraction alveola. The PARASORB® Cone and Fleece HD mould themselves into the site of extraction or excision wound respectively, stops bleeding and will be absorbed.

STORAGE

Store between 0°C - 25°C.

Store protected from moisture and direct sunlight. Do not refrigerate.

Store out of the sight and reach of children.

PRODUCT CONFORMS TO THE ESSENTIAL REQUIREMENTS OF THE COUNCIL DIRECTIVE 93/42/EEC CONCERNING MEDICAL DEVICES.

 **MANUFACTURER**
RESORBA Medical GmbH
Am Flachmoor 16
90475 Nürnberg, Germany
Tel. +49 9128 / 91 15 0

CE1434

GA 0071, Status 2020-05

CONTRAINDICATIONS

The PARASORB® Cone and Fleece HD must not be used in patients with a known hypersensitivity to products of equine origin.

The PARASORB® Cone and Fleece HD must not be used in infected areas.

The PARASORB® Cone and Fleece HD must not be used during pregnancy and lactation.

WARNINGS

Chlorine-releasing antiseptic agents, tannins and caustics should not be used in combination with collagen.

PARASORB® Cone and Fleece HD are not intended to be sutured through; doing so may compromise the performance of the PARASORB® Cone and Fleece HD.

SIDE EFFECTS

Due to the xenogenic origin of the collagen, there is a possible risk of intolerance reactions, despite the confirmed low antigenicity.

FURTHER NOTES

Due to the affinity of collagen to bleeding surfaces, dry instruments and gloves have to be used when applying PARASORB® Cone and Fleece HD.

The product is intended for single use only. PARASORB® Cone and Fleece HD must not be reused, there is a risk of infection for the patient and/or user and the functionality of the product cannot be guaranteed. Remaining material, left over due to cutting for example, must be disposed of according to local procedures and guidelines.

1

2

3

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

-  = REFERENCE NUMBER
-  = BATCH NUMBER
-  = USE BY YEAR – MONTH – DAY
-  = CONSULT INSTRUCTIONS FOR USE
-  = DO NOT RE-USE
-  = DO NOT RESTERILIZE
-  = DO NOT USE IF PACKAGE IS DAMAGED
-  = STERILIZED USING ETHYLENE OXIDE
-  = TEMPERATURE LIMIT
-  = PROTECT FROM DIRECT SUNLIGHT
-  = MANUFACTURER
-  = HIBC-CODE
-  = KEEP DRY

CE1434 = CE MARKING AND IDENTIFICATION NUMBER OF THE NOTIFIED BODY

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