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EZ Cure should be removed from the packaging using sterile gloves or instruments. Mukoperiosteal incision flaps are developed in the site to be treated. The incision should be subcircular as possible. Clinician should perform thorough debridement. As much as possible tissue must be preserved to allow for primary closure of the wound and correct positioning of the flaps.

- Bone defect should be filled without excess with bone grafting material first, MBCP etc..

- EZ Cure can be cut with sterile scissors to desire shape before being hydrated. EZ Cure can be placed either dry or hydrated. For better handling characteristics, the membrane can be hydrated in sterile water or saline solution for approximately not more than five minutes prior to the final placement.

- EZ Cure should overlap the walls of the defect by at least 2mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

EZ Cure is applied over the defect without further treatment and held in place with moderate pressure. The period of time necessary to apply pressure will vary with the degree of bleeding. Adherence to the bone surface is achieved by gel formation of the collagen fibres with blood.

- Complete penetration of the membrane by blood and exudates allows perfect adaptation and adhesion of the membrane to bony structures and makes the formation of a blood clot possible.

- Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination.

- Due to the high tensile strength of this membrane, fixation is possible. Fixation may be indicated to avoid membrane displacement due to loading or mobilization.

- The mucoperiosteal flap is sutured over the collagen membrane (e.g., single sutures and deep mattress sutures).

- The wound should be closed completely to avoid accelerated resorption due to membrane exposure.

Special instructions for use in periodontology:

- A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated.

A postoperative maintenance phase can help to ensure long-term therapeutic success.

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LIMITATIONS FOR USE:

Due to the adherence to the bone tissue and the elasticity of EZ Cure bone augmentation material is required to create and maintain space for bone formation (ex: MBCP...) EZ Cure is therefore not indicated for single use in guided bone regeneration (GBR) without any space-making material. EZ Cure should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated.

The material has not been tested on pregnant women.

POST OPERATIVE CARE

The patient should be monitored closely.

- If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions are recommended.

- In the event that a membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curedtted to remove any inflamed or infected tissue.

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