Hyaloss® matrix

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2011/01

Hyaloss® matrix

HYAFF®-based biopolymer adjunct for bone grafting in dental surgery

Product information

Hyaloss® matrix is a biodegradable gelling adjunct, composed of loose fibers, which has been specifically designed to address handling problems associated with bone grafting.

Hyaloss® matrix is entirely composed of HYAFF®, a derivative of hyaluronic acid, a naturally occurring extracellular matrix molecule.

Giving that bone grafts, both autologous and allogeneic, have been demonstrated to be effective in dental bone reconstruction, their use in parodontology is presently considered the gold standard for the treatment of bone voids and for the restoration of the bony structure where this is missing. The application of bone grafts can be improved by using Hyaloss® matrix because of its rapid gelling properties when in contact with blood or sterile saline. When mixed with bone chips it results in a paste that is easy to handle and apply.

The use of gelling fibers to obtain the amalgam allows the surgeon to prepare a bony paste with adequate consistency, improves handling and retention at the site of application.

Hyaloss® matrix is bioabsorbed in few days and does not affect the natural bone reconstruction process.

Characteristics of the device.

Hyaloss® matrix is a fast gelling adjunct composed of fibers made of HYAFF®, an esterified derivative of hyaluronic acid.

The main characteristics of the fibers are the following:

1. Complete biocompatibility: the polymer used is a derivative of hyaluronic acid (one of the main components of the extracellular matrix). After spontaneous degradation, the polymer releases the parent molecule enriching the transplantation site with hyaluronic acid.
2. Ease of handling and application.

Indications for use

Hyaloss® matrix is a biodegradable material for use in dental surgery as an adjunct for bone graft application in bone voids or defects. Hyaloss® matrix is employed to facilitate handling and application of bone chips without affecting the normal process of osseous regeneration made by the bone graft component.

Contraindications

Individual hypersensitivity to the product

Instructions for use

1. Remove Hyaloss® matrix from the primary packaging and place it on a sterile surgical field
2. Cut the fibers, if needed, using sterile scissors
3. Mix fibers with the bone graft to form a paste. The fibers turn into a gel when in contact with blood. If required, sterile saline can be added to help gel formation.
4. The amount of fibers in each container is designed to be mixed with 0.5 ml of bone graft, but is up to the dental specialist's determination of the ratio between the gelling agent and the bone component.
5. Fill the bone defect with the paste obtained and proceed with the standard surgical technique.

6. Hyaloss® matrix should be dissolved in about one week.

Warning and precautions

Hyaloss® matrix is sterilised by gamma-radiation. It is recommended to use Hyaloss® matrix immediately after the opening of the tray, discard any unused portion of the device. There is a risk of non sterile product and patient infection, in case any unutilized portions of product were to be re-used after first application.

The product is for single use only: sterility is guaranteed as long as the package is closed and undamaged.

Do not re-sterilize the product. FAB assumes no liability for products that have been re-sterilized by health care facilities.

Do not use in case of package damage and report to local distributor.

After use, dispose according to applicable national practice.

Store in a cool and dry place (T < 40°C).

Do not use after the expiry date reported on the box and on the primary packaging.

Keep out of reach of children.

Presentation

4 individually packaged trays each containing fibrous adjunct for 0.5 cc bone graft.

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Anika Therapeutics S.r.l.