
# BetaBASE<sup>®</sup> MP

**BIOVISION**  
BIOMATERIAL

BIOVISION GmbH  
Am Vogelherd 52  
98693 Ilmenau  
Germany  
info@biovision.de  
+49 3677 64 07 0  
+49 3677 64 07 13

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## EN INSTRUCTION MANUAL

### 1. Specific purpose

BetaBASE<sup>®</sup> MP is a bioresorbable bone replacement material for solid tissue replacement. It serves for the filling of mechanically largely non-stressed, pathological, traumatic or surgery-related bone defects in oral and maxillo-facial surgery. Before using it is recommended to verify the suitability of the products for the planned intervention.

#### Indications

- Filling of defects after extirpation of bone cysts
- Filling of extraction defects for creation of an implant bed
- Filling of defects with correction osteotomies
- Other multiple-wall bone defects of the alveolar processes and the facial bones
- Periodontal defects, also in association with membranes
- Defects after extraction of retained teeth
- Defects after apicoectomies
- Sinus floor elevation
- Filling of gaps between alveoli and implants

#### Contraindications

- Medication and therapies that influence bone wound healing, e.g.: Diabetes mellitus, treatments with glucocorticoids and antineoplasm acting agents or irradiation with  $\gamma$  radiation
- Osteomyelitis
- Infections at the site of the implantation

### 2. Factors that can impair the success of the operation

- Insufficient bone quality (e.g. osteoporosis)
- Local bone tumours
- Systemic diseases and metabolic disorders
- Anamnestic indications for infections
- Drug dependency, alcohol and medication abuse
- Mental inability of the patient to understand and follow the instructions given by the physician

#### Possible negative effects

With the application of BetaBASE<sup>®</sup> MP it can in rare cases come to the occurrence of the following product, respectively procedure related negative effects.

- Infections
- Bone fractures or pseudarthrosis of an existing fracture
- Bleeding events
- Delayed/disturbed wound healing
- Injuries that are caused through the implantation, including vascular injuries, nerve lesions and tissue injuries
- Heterotopic ossifications

### 3. Patient briefing

The physician is to inform the patient about possible impairments of the operation success and possible negative effects according to the sections “contraindications” and “possible negative effects”, insofar as applicable for the patient. Furthermore, the patient is to be informed, which measures he or she can undertake to reduce possible impairments and negative effects.

### 4. Operation technique

A surgeon that is not familiar with the specific product is expected to exactly study the surgery technique before the operation.

### 5. Handling

**NOTE:** This product is packed in a protective and sterile packaging.

**WARNING:** Danger of injury and risk of damage of products: The non-observance of these instructions for use can lead to injuries to patients, users and third parties, as well as to damage of the product. Carefully read through the instructions for use and always adhere to the described instructions. Verify the function of the products used in combination.

**WARNING:** Before using verify the expiry date, as well as the packing for any damage. Products with overdue expiry date, respectively damaged packaging must not be used.

**WARNING:** Risk of infection: Through the use of non-sterile implants there is risk of infection for the patient. Only remove the implants from the sterile packing immediately before the application.

**WARNING:** Risk of infection: A reconditioning of the disposable article can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

### 6. Pharmaceutical form / Package size

BetaBASE<sup>®</sup> MP is available as follows:

Article number	Particle size	Package size
BMP 19002105	0,2 – 0,5 mm	1 x 0,5 ml
BMP 19002505	0,2 – 0,5 mm	5 x 0,5 ml
BMP 19005105	0,5 – 1,0 mm	1 x 0,5 ml
BMP 19005110	0,5 – 1,0 mm	1 x 1,0 ml
BMP 19005505	0,5 – 1,0 mm	5 x 0,5 ml
BMP 19005510	0,5 – 1,0 mm	5 x 1,0 ml
BMP 19010105	1,0 – 2,0 mm	1 x 0,5 ml
BMP 19010110	1,0 – 2,0 mm	1 x 1,0 ml
BMP 19010505	1,0 – 2,0 mm	5 x 0,5 ml
BMP 19010510	1,0 – 2,0 mm	5 x 1,0 ml

### 7. Application

BetaBASE<sup>®</sup> MP should be exclusively used by skilled users with corresponding training and experience in implantology and augmentation technique. It is recommended to use the granulate in moist condition, the mass of the individual particles adheres in moist condition and allows more accurate insertion into the defect than in dry form. For moistening sterile isotonic saline solution can be used, as well as own blood of the patient or blood plasma (PRP), with sus-

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## EN INSTRUCTION MANUAL

pected infections under addition of an antibiotic solution. If the defect is larger than 2 cm<sup>3</sup>, we recommend mixing BetaBASE<sup>®</sup> MP with autologous bone material, if the removal is anatomically possible. During the application, the material should not be compressed but evenly layered into the defect, respectively filled to maintain the interconnecting porosity. The selection of the particle size to be used generally depends on the size of the defect. Thereby the following rule of thumb applies: The larger the size of the defect, the larger the selected particle size should be. Depending on the size the introduction into the defect can be performed with spoon or spatula. After the filling of the bone defect, which should take place volume filling up to the edges, the wound is closed saliva-proof with the usual suturing technique. The covering gingiva-periost should achieve a primary closure without tension, so that no granulate particles can be washed out. If a primary closure is not possible, depending on the situation we recommend the augmentation with the non-resorbable Cytoplast membrane under observation of the instructions for use for the membrane technique. During the initial healing phase the augmented area should not be subjected to stress (e.g. through contact with a temporary prosthesis etc.). Depending on the regeneration dynamics of the patient, the resorption time is approximately 6 to 24 months.

### 8. Traceability

The self-adhesive identification labels that are enclosed with the original packing must be stored for traceability and be adhered to the operation report, as well as in the patient file. The individual responsible for the stock-keeping must maintain records of the available products, including article number (REF) and batch number (LOT), so that in the event of changes or recalls action is quickly possible.

### 9. Packaging and sterility

Every individual product is double-packed and sterilised through gamma radiation with a dose of at least 25 kGy. The sterility is guaranteed up to the expiry date, insofar as the packing is neither opened nor damaged.

### 10. Product material

BetaBASE<sup>®</sup> MP is a pure  $\beta$ -tricalcium phosphate ceramic. In the individual granulate particles there are micro pores (up to 5  $\mu$ m), as well as also macro pores (up to 1 mm).

Composition:

Calcium oxide CaO 52.0 to 54.2 weight %  
Phosphorus pentoxide P<sub>2</sub>O<sub>5</sub> 45.8 to 48.0 weight %  
The ratio Ca/P is 1.5.

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### 12. Storage conditions

The product must be stored dry in the unopened original packing. Products with elapsed expiry date, respectively damaged packing must not be used.

### 13. Responsibility

Complications or other negative effects due to incorrect indication or surgical technique, unsuitable selection or handling of products or instruments, asepsis, etc. fall under the responsibility of the surgeon.

### 14. Conformity with directives

According to the medical device directive (MDD) 93/42/EEC this medical device is labelled with the CE label. If there is an identification number after the CE label then this designates the responsible specified authority.

### 15. Disposal

No special measures are necessary during the disposal of this medical device. The national laws and directives must be observed.

### 16. Packing symbols

	Manufacturer		Do not reuse
	Observe instructions for use		Do not re-sterilise
	Use before		Batch number
	Do not use product with damaged packing		Article number
	Sterilised through irradiation		Store dry

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