Instruction for Use (ENG)

DESCRIPTION

In^{oss}™ is a Resorbable Bone Graft Substitute that is rapidly replaced by newly-formed bone. This material is a non self setting providing mechanical resistance after healing and osseous rehabilitation. In^{oss}™ is composed of biphasic Calcium Phosphate (60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP)) particles associated with an excipient vehicle of pharmaceutical grade quality (hydroxypropylmethylcellulose in an aqueous solution). The soluble carrier acts as a spacer and a binder of the particles. In^{oss}™ provides an environment for bone ingrowth. On time In^{oss}™ is fully resorbable. In^{oss}™ is provided sterile for single patient use in a syringe.

INDICATIONS

In^{oss}™ is intended for use as a bone grafting material to fill, augment or reconstruct osseous bone defects in particular in periodontal or oral/maxillofacial applications. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. In^{oss}™ can generally be used for bone filling in closed cavities. In^{oss}™ can be used with autogenous bone grafting materials.

Typical uses include but are not limited to:

- Periodontal/Infrabony defects
- Ridge augmentation
- Extraction sites (implant preparation/placement)
- Sinus lifts
- Cystic cavities

CONTRA-INDICATIONS

In^{oss}™ has limited initial mechanical properties. Therefore, this product is contraindicated where the device is intended as structural support. Conditions representing contraindications include also:

- Degenerative bone disease
- Infectious diseases (evolution or uncontrolled)
- Immune deficiency (congenital or acquired)
- As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure. This includes (but is not limited to) individuals with long-term steroid therapy or treatment acting on the calcium or phosphorus metabolism.
- Oslerian risk
- Unbalanced diabetes
- Opening mening
- Cellulitis

- Local:
 - Trophic trouble
 - Irradiated patient
 - Osteomyelitis acute or chronic
 - Alveolitis dry or infectious not treated
 - Granuloma or cyst not cured
 - Intra-articular implantation
 - Opening mening
 - Cellulitis
 - Implantation in necrotic surgical sites

ADVERSE-EFFECTS

Possible adverse effects include but are not limited to:

- Wound complications including haematoma, infection, and other complications that are possible with any surgery

WARNING AND PRECAUTIONS

- In^{oss}™ is intended only for bony voids or gaps that are not intrinsic to the stability of the bony structure.
- Maximum contact between the product and the recipient bone must be established.
- Do not overfill defects.
- In^{oss}™ does not harden after implantation until bone rehabilitation. The device should be secured to prevent motion and migration, use in areas where the graft can be adequately contained.
- Do not leave defect open, (e.g. mucosa well sutured and/or use of appropriate membrane).
- Do not compromise blood supply to the defect area.
- In^{oss}™ has to be in direct contact of at least 2 or 3 surfaces. Maximum contact between In^{oss}™ and recipient bone must be achieved.
- In^{oss}™ does not harden after implantation until healing. In^{oss}™ is a void filler used in cavity with respected corticals or associated with an osteosynthesis device. To avoid any possibility of leakage out of the cavity of implantation, bone flap or resorbable membrane can be associated.
- Adequate post-operative immobilization is to be used with In^{oss}™ to prevent movement that may lead to soft tissue ingrowth.
- Until healing is complete In^{oss}™ may not withstand weight bearing. Post surgery discharge is necessary.
- Implantation is not recommended in revision surgery site containing non resorbable biomaterials residues. Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant are important considerations in the successful utilization of this device.
- The product has not been tested on pregnant women. The risk of health has not been established.

STERILIZATION AND CONSERVATION

In^{oss}™ is supplied sterile in a syringe. In^{oss}™ is sterilized by steam. This product should not be resterilized. In all operations involving In^{oss}™, appropriate aseptic technique must be observed. In^{oss}™ should be stored dry at ambient temperature and should not exceed 50°C (122°F). This device is for single patient use only and should never be resterilized. The sterilization process of this product was validated during its development. We do not guarantee the sterility of the product if a new sterilization is performed on it. The non-respect of the single use can conduct to septic issues. Additionally, the resorption process of this bone substitute started immediately after its implantation on the osseous defect so it can not be re-used.

DO NOT USE IF PACKAGE IS OPENED OR DAMAGED OR IF EXPIRATION DATE HAS BEEN EXCEEDED.

This product is to be handled or implanted by trained qualified persons having read this instruction for use.

Caution: Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

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